



With over 100 years of commitment to quality, we prioritise safety and reliability. To meet the diverse needs of our customers, we offer a comprehensive range of products.

# Daichan Single Jaw type



The Daichan series allows you to complete working models by attaching impression models to the base (Daichan) with dowel pins and simply separating them. Compared to conventional methods, it reduces the time to create detachable models to one-third, significantly saving time and cost.

- Time-saving No need to plant dowel pins or create retention grooves.
- High Precision Produces dense models unaffected by the setting expansion of gypsum.
- Eco-friendly No need to mix gypsum, reducing industrial waste.
- CAD/CAM Compatible No plastic parts reduce scanning errors due to reflections.

## Daichan



50 pieces per pack Gray × Yellow

### Features

- Equipped with 5 dowel pins arranged at equal intervals.
- Suitable for single crowns.

## Daichan II



50 pieces per pack Gray × White

### Features

- Equipped with 6 dowel pins arranged at staggered intervals.
- It is suitable not only for single crowns but also for multiple crowns.
- It is equipped with wings for mounting on the articulator.

### <Daichan, Daichan II Usage instructions>

- ①Smooth the base of the dental model with a trimmer. ※Ensure the base is flat to avoid detachment during adhesion.
- ②Mark the areas where you will insert the saw (lines).
- ③Check the position of the “abutment teeth” and the pins before setup.  
If the pins are not visible, use the anti-rotation groove as a guide.
- ④Apply a thin, even layer of gel-type instant adhesive (Daichan Bond) to the base of the model.
- ⑤Attach the model to the Daichan, aligning the pin positions.
- ⑥Solidify any excess adhesive with an accelerator.
- ⑦Confirm the adhesive has cured and separate the model using a saw to complete.

## Related Products

## Daichan Bond



20g per bottle

Adhesive for model base

### Features

- Features strong adhesion and ample curing time.
- Gel-like consistency prevents dripping, making it usable on vertical surfaces.

For inquiries and product orders, please contact the following

Visit our website!



Manufacturer and Distributor


**SAN ESU GYPSUM CO.,LTD**

<https://www.sanesugypsum.com/dental/>



With over 100 years of commitment to quality, we prioritise safety and reliability. To meet the diverse needs of our customers, we offer a comprehensive range of products.

# Daichan Full Jaw type



The Daichan series allows you to complete working models by attaching impression models to the base (Daichan) with dowel pins and simply separating them. Compared to conventional methods, it reduces the time to create detachable models to one-third, significantly saving time and cost.

- **Time-saving** No need to plant dowel pins or create retention grooves.
- **High Precision** Produces dense models unaffected by the setting expansion of gypsum.
- **Eco-friendly** No need to mix gypsum, reducing industrial waste.
- **CAD/CAM Compatible** No plastic parts reduce scanning errors due to reflections.

## / Daichan Full Jaw



### Features

- Equipped with 14 dowel pins.
- Uses single-sided cut pins.

30 pieces per pack Gray × Yellow

### <Daichan Full Jaw, Daichan Full Jaw L Usage instructions>

- ① Smooth the base of the dental model with a trimmer.  
※ Ensure the base is flat to avoid detachment during adhesion.
- ② Mark the areas where you will insert the saw (lines).
- ③ Check the position of the “abutment teeth” and the pins before setup.  
If the pins are not visible, use the anti-rotation groove as a guide.
- ④ Apply a thin, even layer of gel-type instant adhesive (Daichan Bond) to the base of the model.
- ⑤ Attach the model to the Daichan, aligning the pin positions.
- ⑥ Solidify any excess adhesive with an accelerator.
- ⑦ Confirm the adhesive has cured and separate the model using a saw to complete.

- **Easy Separation** Uses a retention method without dowel pins.
- **Low Thermal Deformation** Uses a special resin with heat resistance up to 250°C for the plastic plate.

## / Pla Daichan F




30 pieces per pack Gray × Yellow

### <Pla Daichan F Usage instructions>

- ① Smooth the base of the dental model with a trimmer. ※ Ensure the base is flat to avoid detachment during adhesion.
- ② Confirm the base is flat and trim the model to fit the Pla Daichan F.
- ③ Apply a thin, even layer of gel-type instant adhesive (Daichan Bond) to the base of the model.
- ④ Attach the model to the Pla Daichan F, aligning the pin positions.
- ⑤ Solidify the excess adhesive with an accelerator as soon as it overflows.
- ⑥ Verify the curing of the adhesive, and then separate the model using a saw.
- ⑦ Gently tap the base with a small hammer to remove the dental base.

Manufacturer and Distributor  
**SAN ESU GYPSUM CO.,LTD**



<https://www.sanesugypsum.com/dental/>


## / Daichan Full Jaw L type



### Features

- Equipped with 14 dowel pins.
- Expanded plate size increases applicable cases.
- Uses double-sided cut pins for improved anti-rotation effect.

30 pieces per pack Gray × White



For inquiries and product orders, please contact the following