


# Other dental-related products

With over 100 years of commitment to quality, we prioritise safety and reliability. To meet the diverse needs of our customers, we offer a comprehensive range of products.



## Phosphate-Based Dental Investment

### My Vest · My Vest Special Liquid

Common name: My Vest High-Temperature Dental Casting Investment

Notification number: 28B2X1003700033

Common name: My Vest Special Liquid Auxiliary Material for High-Temperature Dental Models

Notification number:

My Vest is a phosphate-based high-temperature dental casting investment developed for casting high-melting-point alloys. It is suitable for casting non-precious alloy crowns, bridges, inlays, and metal bonds.



My Vest 3kg Aluminum Bag  
My Vest Special Liquid 500ml

#### Example of Special Liquid Dilution

Ring Casting	Mix for 30 seconds	leave at room temperature for 40 minutes	Special Liquid 15ml+Water 3ml
Ringless Casting	Mix for 30 seconds	leave at room temperature for 40 minutes.	Special Liquid 13ml+Water 5ml
Ring Casting	Mix for 1 minute	leave at room temperature for 1 hour	Special Liquid 12ml+Water 6ml

※The most standard usage method is used as a basis. Please make slight adjustments to the optimal conditions according to the usage conditions.

#### Fabrication of non-precious alloy crowns, bridges, inlays, and metal bonds

Liquid-to-Powder Ratio	0.18
Setting Time	12 min
Linear Setting Expansion	0.60 %
Compressive Strength	4.8 MPa

(Measured under the conditions of room temperature 22°C, water temperature 22°C, 100g of powder, 15ml of the dedicated solution, and 3ml of water.)

#### Features

- Uses a system that adjusts the concentration of the special liquid according to the casting method, providing excellent adaptability.
- Ensures precise castings with high permeability.
- Facilitates easy removal after casting.
- Produces a smooth casting surface, reducing the time required for polishing.

## CAD/CAM-related dental products

### Sun Ace Disk W

Common name : Dental casting wax

Notification Numbers : 28B2X1003700055

Sun Ace Disk W is a wax disk specialized for CAD/CAM milling, excelling in cutting performance and leaving minimal residue after incineration. It is suitable for rapid heating.



One disk per box

Diameter	9 8 mm	9 8 mm
Thickness	1 4 mm	2 0 mm

#### Wax disk specifically for CAD/CAM milling

#### Features

- High Cutting Performance: Provides a clean finish.
- Minimal Residue After Incineration: Leaves almost no residue.
- Easy-to-See Shading: Comes in a gray color that is easy on the eyes.
- Two Size Options: Available in diameters of 98mm and thicknesses of 14mm and 20mm.
- Additional Build-Up: Allows for additional build-up with margin waxes, etc.

For inquiries and product orders, please contact the following

Visit our website!



Manufacturer and Distributor

**SAN ESU GYPSUM CO.,LTD**

<https://www.sanesugypsum.com/dental/>