



With over 100 years of commitment to quality, we prioritise safety and reliability. To meet the diverse needs of our customers, we offer a comprehensive range of products.

Alginate impression materials • Related Product



Alginate impression materials

New Algistar

Common name: Dental alginate impression materials Certification Numbers: Fast 219AABZX00188000, Normal 219AABZX00189000

Dental alginate impression material available in powder form, with two types differing in working time. Minimizes powder dispersion, easy to mix, with a smooth finish that adapts to various plasters. It does not sag in the oral cavity and is easy to remove after setting, making it friendly for both users and patients.



5kg box (1kg aluminum bags × 5)



Green



5kg box (1kg aluminum bags × 5)




Pink



※Please use a dedicated measuring instrument (sold separately) for weighing.

For taking oral impressions.

Type	Fast	Normal
Working Time	1'30"	2'00"
Intraoral Setting Time	1'30"	1'30"



Features

- **Quick Integration with Water** Sinks quickly in water, making it easy to mix.
- **Easy and Fast Mixing** Can be mixed quickly with light effort, reducing strain on the wrist.
- **Non-Sagging Soft Paste** Does not sag easily due to its sticky paste, even without intentionally making it firm.
- **Almost No Scent** Lacks a strong smell, ensuring it can be used comfortably by patients sensitive to odors.
- **Easy to Remove from Teeth** The set impression material does not adhere strongly to the teeth, making it easy to remove from the oral cavity, reducing patient discomfort.
- **Minimal Deformation** The impression material is flexible, reducing deformation even when removed from areas with significant undercuts.
- **Compatible with Automatic Mixers** Suitable for both hand mixing and automatic mixers.

Tray cleaning material

AL Remover

AL Remover is a liquid tray cleaner that shrinks alginate impression material, allowing it to be easily peeled off. It does not dissolve the alginate impression material, enabling fast processing. It leaves no residue, has a long lifespan, and also offers antibacterial properties.



2000 ml plastic bottle

liquid type

Features

- **Fast** It removes the alginate impression material by shrinking it without dissolving it, allowing for quick removal.
- **Easy** It can be used as is and lasts long without getting dirty.
- **Clean** Since it does not dissolve the alginate impression material, no residue is left, and it has antibacterial properties.

San Esu Tray Cleaner

San Esu Tray Cleaner is a powerful tray cleaner developed with the goal of being more efficient and providing a cleaner result.



1kg poly bottle (comes with one measuring spoon)

powder form type

Features

- **Fast** Quickly dissolves and softens alginate impression materials attached to trays, allowing for efficient cleaning.
- **Easy** Changes color from blue to clear as its effectiveness decreases, making it easy to know when to replace the solution.



※The solution is marine blue.

For inquiries and product orders, please contact the following

Visit our website!



Manufacturer and Distributor

SAN ESU GYPSUM CO.,LTD

<https://www.sanesugypsum.com/dental/>